

Precellys® 24

LYSIS & HOMOGENIZATION AUTOMATED EQUIPMENT



The perfect solution for biological sample preparation to extract nucleic acids (DNA, RNA), proteins

- Up to 24 tubes simultaneously
- Ergonomic loading and cleaning
- Reliability

Rapid and reliable biological sample lysis

Precellys® 24 is a benchtop equipment dedicated to the grinding, lysis and homogenization of biological samples. It can simultaneously operate up to 24 tubes containing up to 2 ml sample at very high speed.

Flexible and easy setting protocol : buffer and samples are added in 2 ml tubes pre filled with specific beads. Times, speed and number of cycles are defined before run.

Ergonomic loading and cleaning : fast and easy loading thanks to the innovative patented *Press'n Block* system. Easy cleaning system with a direct access to the sample holder and compatibility to usual cleaning solutions.

Efficiency : Precellys® 24 is a high efficiency solution for grinding difficult samples like spores of bacteria and hard tissues. The high speed and specific motion guarantee an homogeneous and efficient grinding.

Reliability : specifically developed and designed to be used in L3 laboratory under high throughput conditions, 50 to 100 runs per day.


Cross contamination free



Up to 24 tubes



Patented Press'n Block



Precellys[®] 24 : Lysis and homogenization automated equipment



Ergonomic loading and settings

- Capacity : up to 24 samples per run
- Same kinetic for each tube
- Cross contamination free - single use tube
- Standard tubes of 2ml with ceramic or metal grinding beads references on demand


Display and control panel

- Run time : 5 - 90 s by 5 s increments
- Rotation speed : 5000 - 7000 rpm by 100 rpm increments
- Number of cycles : 1- 3
- LCD with back light up to 2 lines - for button keyboard



Technical Specifications


- Patented Automated System for blocking tubes and motion
- FCC Part 15 and CE compliance
- Developed to be used in L3 Laboratory
- Dimensions : 290 (L) x 375 (W) x370 (H) mm
- 600 mm (H open cover)
- Weight : 24 kg
- Universal Power : 100 - 230 VAC / 50 - 60 Hz



Animal tissues



Vegetal tissues



Micro-organisms





precellys 24

lysis & homogenization

BERTIN TECHNOLOGIES

Headquarter

Parc d'activités du Pas du Lac
10 bis, avenue Ampère - B.P. 284
78180 Montigny-le-Bretonneux
FRANCE

BERTIN TECHNOLOGIES

Post Address

B.P. 284
78053 Saint-Quentin-en-Yvelines Cedex
FRANCE

Phone : +33 (0) 139 306 169 / 070

Fax : +33 (0) 139 306 185

WebSite : www.bertin.fr

Contact : Biotech Systems Department

E-mail : precellys@bertin.fr

Certifié ISO 9001 - 2009 par



Bertin Technologies
Montigny, Aix, Tarnos

D8000.421.SL003.A
2004, november