

General Purpose Media

UltraCULTURE™ Serum-free Medium

- Delivers excellent growth and expansion of a broad range of adherent and suspension cells
- Simplify your media needs, no FBS needed
- The formulation for UltraCULTURE™ has been submitted to the FDA as a product Master File

PC-1™ Chemically Defined, Serum-free Medium

- Superior media for primary adherent cells

HL-1™ Chemically Defined, Serum-free Medium

- Unique media for hybridoma and lymphocytic cells

	General use	Non-animal origin	Protein-free	Chemically defined	Optimized for
UltraCULTURE™ medium	■	–	–	–	Most cell types
PC-1™ CD medium	■	–	–	■	Adherent primary cells and cell lines
HL-1™ CD medium	–	–	–	■	Hybridomas and many adherent cell lines

Ordering information

Cat. no.	Product	Size
12-725F	UltraCULTURE™ medium	500 ml
77232	PC-1™ CD medium	2 X 500 ml
77201	HL-1™ CD medium	2 X 500 ml